




Inspecting Eyewear using the ALM-700/800



1. Set the lens meter to the complete pair single vision setting:

1A: Use the Measurement Selection icon to change lens type.  


1B: Press the Clear icon to erase previous readings. 

1C: Press the Lens icon to change from single lens to complete pair.  



2. Place the right eye over the lens stop:

2A: Position the bottom of the frame so its resting against the frame table with the right-side nose pad/bridge in contact with the left nose pad that's attached to the frame table.



3. Move the cross to the center of the target:

3A: While keeping the frame against the frame table, center the cross within the target. Use your left hand to move the cross right or left and the frame table lever to move it up or down.


3B: Make fine movements to get the cross in the exact center of the target. As the cross approaches the center, the screen will display Alignment OK. When the cross is centered it will change color from pink to blue and the screen will display Marking OK.



4. Clamp the lens:


4A: After the cross is centered in the target and the screen displays Marking OK, clamp the R lens in place.

4B: When the lens is clamped, the cross may move a little out of alignment and need to be readjusted. With one hand, release just a bit of pressure from the lens clamp, and readjust the alignment cross so it's perfectly centered and re-clamp.




5. Press the Memory Button:

5A: When the lens is perfectly placed and clamped, press the Memory Button. The color of the measurement value box will change, indicating the measurement values are saved.




6. Mark the lens:

6A: After the measurements are saved, use the Marking Pin Lever to spot the optical center of the lens. This marking will be used later in the process to inspect OC Height.




7. Remove the right lens:

7A: Lift the lens clamp and remove the R lens.



8. Repeat Steps 2 – 7 for the left lens:

8A: After the measurement of R lens, place the L lens on the lens stop with the left side nose pad/bridge in contact with the right nose pad that's attached to the frame table, then hold the lens in place with the lens clamp. At this time, the measurement box automatically switches from R lens to L lens based on the position of the fixed nose pad.



9. Compare Sphere power, Cylinder power, & Axis readings to the Quality Standards:

9A: Compare the recorded values for Sphere power, Cylinder power, & Axis to the allowable tolerances documented on the RxO work ticket or Quality Standards.



Note: Use the "Inspection Tolerance" for optimized prescriptions from RxO.


INSP. TOLERANCE				
Sphere	Tol	Cyl	Tol	Axis
R: +1.75	0.13	-2.00	0.13	180
L: +2.00	0.13	-2.00	0.13	170

Optical Standards - US & Canada		EYEWEAR INSPECTION JOB AID	
Standard	Requirement	Inspection	Pass/Fail
1	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
2	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
3	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
4	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
5	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
6	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
7	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
8	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
9	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
10	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
11	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
12	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
13	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
14	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
15	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
16	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
17	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
18	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
19	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail
20	Prescription must be within tolerance of the RxO work ticket or Quality Standards.	Compare recorded values to tolerance.	Pass/Fail

10. Measure the PD Using the EIP JobAid:

10A: Use the Alignment Grid on the EIP Job aid by centering the bridge and reading off the scale.


Note: the image to the right is an example of a Right eye with the white OC mark reading 30 mm on the scale.



10B: Determine the values on the Right and Left Lenses.

10C: Compare values to the work ticket- If the values match, continue to checking OC Height


If the values do not match continue to the next step



10-1. Determine the power of the lenses at 180°:

10-1A: Use the Lab Ticket to determine the power of the lens at 180- The image on the right shows the Right lens is over 2.75. Go to the next step.

10-1B: If the power is less than 2.75 go to Checking Lenses under 2.75




10-2a. Checking Lenses over 2.75D

10-2A: Compare the value listed on the Lab Ticket to the value you read using the EIP AlignmentGrid:

10-2B: Is the Right OC within 1.5 mm of the prescribed value? If YES, the lens passes.

10-2C: Is the Left lens within 1.5 mm? If YES, the lens passes.


10-2D: Is the pair within a total of 2.5 mm from prescribed? If YES, the eyewear passes inspection for PD Continue to the Check OC Height.



10-2b. Checking PD in lenses under 2.75 D:

10-2A: 1. Use an AR pen to mark where the OC should be both horizontally and vertically.

For example- the lens below marked in the lensometer at 30 mm. The lab ticket listed 32 as the correct PD. Mark a cross at 32 mm




10-2B: Center the cross in the Lens Stop . Use the Axis Marking Pins to ensure the middle pin is directly over the center of the cross.

10-2C: Clamp the lens and press the memory button to capture the prism value.

10-2D: The horizontal prism values may be up to 0.33D in each eye.


10-2E: If both lenses are equal to or less than 0.33, for a total of 0.67, the pair passes inspection for PD



In tolerance
below .33D

11. The OC's were spotted in Step 1. If needed re-spot by:

- 11A: Align the Right Lens until the Target is centered and the screen reads "Marking Ok"
- 11B: Use the marking pins to mark the Right OC
- 11C: Align the Left lens until the Target is centered and the screen reads "Marking OK"
- 11D: Use the marking pins to mark the Left OC



12. Measure the OC Height Using the EIP JobAid:


12A: Use the Alignment Grid on the EIP Job aid by aligning the frame to the dark bottom line and reading off the scale.

Note: the image to the right is an example of a Right eye with the OC mark reading 20 mm on the scale.

12B: Determine the values on the Right and Left Lenses.

12C: Compare to the work ticket- If the values match the pair passes.

If the values do not match continue to the next step



12-1. Determine the power of the lenses at 90°:

12-1A: Use the Lab Ticket to determine the power of the lens at 90- The image on the right shows the Right lens is over 3.75. Go to the next step.

12-1B: If the power is less than 3.75 go to Checking Lenses under 3.75

Power @	
90	180
- 4.42	- 4.33
- 4.99	- 4.51
Vendor	

12-2a. Checking Lenses over 3.75D:

12-2A: Compare the value listed on the Lab Ticket to the value you read using the EIP Alignment Grid:

12-2B: Is the Right OC within 1.0 mm of the prescribed value? If YES, the lens passes.


12-2C: Is the Left lens within 1.0 mm? If YES, the lens passes.

12-2D: Are they within 1.0 mm of each other? If YES, the pair passes inspection for OC Height

12-2b. Checking OC's in lenses under 3.75 D:

12-2A: Use an AR pen to mark where the OC should be both horizontally and vertically.

For example- the lens below marked in the lensometer at 18 mm. The lab ticket listed 20 as the correct OC. Mark a cross at 20 mm.



12-2B: Center the cross in the Lens Meter. Use the Axis Marking pins to ensure the middle pin is directly over the center of the cross.

12-2C: Clamp the lens and press the memory button to capture the prism value.

12-2D: The values may be up to 0.33D in each eye. 17-2E: If both lenses have prism check the direction. If both down or both up and each lens passed individually the pair passes. If one lens has UP prism and one DOWN- add the values together, the result must be equal to or less than 0.33D.




In tolerance
below .33D

1. Recreate the manufacturers lens markings using the EIP job aid:

1A: The following markings are required to final inspect progressive lenses; DRP, NRP, Fitting Cross, PRP and 180-degree engravings.

1B: Determine the lens type from the Job Ticket and recreate the progressive markings using the appropriate layout from the EIP Job aid.



2. Set the lens meter to Progressive lens mode:

2A: Use the Measurement Selection icon to change lens type to progressive.




2B: Press the Lens icon to change from single lens to complete pair.



3. Center the Distance Reference Point (DRP) over the lens stop:

3A: Position the bottom of the frame so its resting against the frame table with the nose pad/bridge in contact with the left nose pad that's attached to the frame table, then center the lens stop within the DRP on the R lens.


3B: Secure the R lens in place with the lens clamp.



4. Press the Memory button to save the measurements:

4A: Make sure the R lens is clamped and the frame is resting against the table before measuring it, unclamped lenses will read inaccurately.

4B: A measurement can be taken even if the messages of Alignment OK and Marking OK are not displayed.




5. Unclamp and position the lens over the Near Reference Point (NRP):

5A: Release the lens clamp and use the frame table lever to move the table towards yourself until the NRP is over the lens stop, then re-clamp the lens.


5B: Use your right hand to move the frame table lever and your left hand to keep the lens in position against the lens stop and the bottom of the frame against the frame table.

5C: If some of the NRP has been cut off and cannot be centered over the lens stop, position as much of the NRP as possible over the lens stop.




6. Press the Memory button to capture the Add power:

6A: Make sure the lens is clamped and the frame is resting against the table before capturing the Add power, unclamped lenses will read inaccurately.




7. Remove the right lens:

7A: Lift the lens clamp and remove the R lens.



8. Repeat Steps 3 – 7 for the left lens:

8A: After the measurement of R lens, position the L lens DRP over the lens stop with the left side nose pad/bridge in contact with the right nose pad that's attached to the frame table. At this time, the measurement box automatically switches from R lens to L lens based on the position of the fixed nose pad.



11-1D: Is the pair within a total of 2 mm from prescribed? If YES, the eyewear passes inspection for PD. Continue to the Check OC Height.

12. Measure the Fitting Cross Height Using the EIP JobAid:

12A: Use the Alignment Grid on the EIP Job aid by aligning the frame to the dark bottom line and reading off the scale.

Note: the image to the right is an example of a Right eye with the Seg Ht. mark reading 20 mm on the scale.

12B: Determine the values on the Right and Left Lenses.

12C: Compare to the work ticket- If the values Match, continue to Prism inspection.

If the values do not match continue to the next step



12-1. Determine the allowable Fitting Cross tolerance.

12-1A: Compare the value listed on the Lab Ticket to the value you read using the EIP Alignment Grid:

12-1B: Is the Right lens within 1mm of specified Fitting Cross height? If YES, the lens passes.

12-1C: Is the Left lens within 1 mm? If YES, the lens passes.

12-1D: Is the pair within a total of 1 mm difference of prescribed? If YES, the eyewear passes inspection for Fitting Cross. Continue to inspecting Prism.

Optical Standards - US & Canada		EYEWEAR INSPECTION JOB AID																																																																																									
<p>Table 1: Fitting Cross Height Tolerance</p> <table border="1"> <thead> <tr> <th>Prescription</th> <th>Right Eye</th> <th>Left Eye</th> <th>Total Difference</th> </tr> </thead> <tbody> <tr> <td>SPH < +1.00</td> <td>±1.00</td> <td>±1.00</td> <td>±2.00</td> </tr> <tr> <td>+1.00 ≤ SPH < +2.00</td> <td>±0.75</td> <td>±0.75</td> <td>±1.50</td> </tr> <tr> <td>+2.00 ≤ SPH < +3.00</td> <td>±0.50</td> <td>±0.50</td> <td>±1.00</td> </tr> <tr> <td>+3.00 ≤ SPH < +4.00</td> <td>±0.25</td> <td>±0.25</td> <td>±0.50</td> </tr> <tr> <td>+4.00 ≤ SPH < +5.00</td> <td>±0.12</td> <td>±0.12</td> <td>±0.25</td> </tr> <tr> <td>+5.00 ≤ SPH < +6.00</td> <td>±0.06</td> <td>±0.06</td> <td>±0.12</td> </tr> <tr> <td>+6.00 ≤ SPH < +7.00</td> <td>±0.03</td> <td>±0.03</td> <td>±0.06</td> </tr> <tr> <td>+7.00 ≤ SPH < +8.00</td> <td>±0.01</td> <td>±0.01</td> <td>±0.02</td> </tr> <tr> <td>+8.00 ≤ SPH < +9.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> <tr> <td>+9.00 ≤ SPH < +10.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> </tbody> </table>		Prescription	Right Eye	Left Eye	Total Difference	SPH < +1.00	±1.00	±1.00	±2.00	+1.00 ≤ SPH < +2.00	±0.75	±0.75	±1.50	+2.00 ≤ SPH < +3.00	±0.50	±0.50	±1.00	+3.00 ≤ SPH < +4.00	±0.25	±0.25	±0.50	+4.00 ≤ SPH < +5.00	±0.12	±0.12	±0.25	+5.00 ≤ SPH < +6.00	±0.06	±0.06	±0.12	+6.00 ≤ SPH < +7.00	±0.03	±0.03	±0.06	+7.00 ≤ SPH < +8.00	±0.01	±0.01	±0.02	+8.00 ≤ SPH < +9.00	±0.00	±0.00	±0.00	+9.00 ≤ SPH < +10.00	±0.00	±0.00	±0.00	<p>Table 2: Prism Tolerance</p> <table border="1"> <thead> <tr> <th>Prescription</th> <th>Right Eye</th> <th>Left Eye</th> <th>Total Difference</th> </tr> </thead> <tbody> <tr> <td>SPH < +1.00</td> <td>±0.25</td> <td>±0.25</td> <td>±0.50</td> </tr> <tr> <td>+1.00 ≤ SPH < +2.00</td> <td>±0.12</td> <td>±0.12</td> <td>±0.25</td> </tr> <tr> <td>+2.00 ≤ SPH < +3.00</td> <td>±0.06</td> <td>±0.06</td> <td>±0.12</td> </tr> <tr> <td>+3.00 ≤ SPH < +4.00</td> <td>±0.03</td> <td>±0.03</td> <td>±0.06</td> </tr> <tr> <td>+4.00 ≤ SPH < +5.00</td> <td>±0.01</td> <td>±0.01</td> <td>±0.02</td> </tr> <tr> <td>+5.00 ≤ SPH < +6.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> <tr> <td>+6.00 ≤ SPH < +7.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> <tr> <td>+7.00 ≤ SPH < +8.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> <tr> <td>+8.00 ≤ SPH < +9.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> <tr> <td>+9.00 ≤ SPH < +10.00</td> <td>±0.00</td> <td>±0.00</td> <td>±0.00</td> </tr> </tbody> </table>		Prescription	Right Eye	Left Eye	Total Difference	SPH < +1.00	±0.25	±0.25	±0.50	+1.00 ≤ SPH < +2.00	±0.12	±0.12	±0.25	+2.00 ≤ SPH < +3.00	±0.06	±0.06	±0.12	+3.00 ≤ SPH < +4.00	±0.03	±0.03	±0.06	+4.00 ≤ SPH < +5.00	±0.01	±0.01	±0.02	+5.00 ≤ SPH < +6.00	±0.00	±0.00	±0.00	+6.00 ≤ SPH < +7.00	±0.00	±0.00	±0.00	+7.00 ≤ SPH < +8.00	±0.00	±0.00	±0.00	+8.00 ≤ SPH < +9.00	±0.00	±0.00	±0.00	+9.00 ≤ SPH < +10.00	±0.00	±0.00	±0.00
Prescription	Right Eye	Left Eye	Total Difference																																																																																								
SPH < +1.00	±1.00	±1.00	±2.00																																																																																								
+1.00 ≤ SPH < +2.00	±0.75	±0.75	±1.50																																																																																								
+2.00 ≤ SPH < +3.00	±0.50	±0.50	±1.00																																																																																								
+3.00 ≤ SPH < +4.00	±0.25	±0.25	±0.50																																																																																								
+4.00 ≤ SPH < +5.00	±0.12	±0.12	±0.25																																																																																								
+5.00 ≤ SPH < +6.00	±0.06	±0.06	±0.12																																																																																								
+6.00 ≤ SPH < +7.00	±0.03	±0.03	±0.06																																																																																								
+7.00 ≤ SPH < +8.00	±0.01	±0.01	±0.02																																																																																								
+8.00 ≤ SPH < +9.00	±0.00	±0.00	±0.00																																																																																								
+9.00 ≤ SPH < +10.00	±0.00	±0.00	±0.00																																																																																								
Prescription	Right Eye	Left Eye	Total Difference																																																																																								
SPH < +1.00	±0.25	±0.25	±0.50																																																																																								
+1.00 ≤ SPH < +2.00	±0.12	±0.12	±0.25																																																																																								
+2.00 ≤ SPH < +3.00	±0.06	±0.06	±0.12																																																																																								
+3.00 ≤ SPH < +4.00	±0.03	±0.03	±0.06																																																																																								
+4.00 ≤ SPH < +5.00	±0.01	±0.01	±0.02																																																																																								
+5.00 ≤ SPH < +6.00	±0.00	±0.00	±0.00																																																																																								
+6.00 ≤ SPH < +7.00	±0.00	±0.00	±0.00																																																																																								
+7.00 ≤ SPH < +8.00	±0.00	±0.00	±0.00																																																																																								
+8.00 ≤ SPH < +9.00	±0.00	±0.00	±0.00																																																																																								
+9.00 ≤ SPH < +10.00	±0.00	±0.00	±0.00																																																																																								

13. Inspect Progressive Prism:

13A: Press the Clear icon to erase previous readings.


13B: The manufacturer's markings, specifically the Prism Reference Point (PRP) is required for this step. If needed, you can remark the lens using the EIP job aid.

13C: Use the Measurement selection icon to change lens type back to SV.




14. Center the right lens Prism Reference Point (PRP) over the lens stop:

14A: Center the PRP on the R lens over the lens stop and clamp the lens in place.




15. Confirm the lens is properly placed at the PRP:

15A: With the bottom of the frame resting against the frame table, verify the R lens PRP is properly placed by using the Marking Lever to spot the lens. The center dot must be superimposed over the PRP for accurate Prism readings.




16. Press the Memory button to capture the Prism readings:

16A: Make sure the R lens is clamped and the frame is resting against the table before pressing the Memory button, unclamped lenses will read inaccurately.




17. Unclamp and remove the right lens:


17A: Release the lens clamp and remove the R lens.



18. Repeat Steps 10 – 16 for the left lens:

18A: Press  then  icons to switch the lens meter to read the left lens.

18B: Repeat steps 10 – 16 for the L lens.



19. Compare the horizontal and vertical prism readings to the Quality standards:

9-14A: Compare the In & Out prism readings to the allowable tolerance for Horizontal Prism listed in the Quality Standards on the EIP job aid.

9-14B: Compare the Up & Down prism readings to the allowable tolerance for Vertical Imbalance listed in the Quality Standards on the EIP job aid.

Optical Standards - US & Canada		EYEWEAR INSPECTION JOB AID	
<p>Horizontal Prism</p> <p>1. Horizontal Prism: The sum of the horizontal prism readings (In & Out) for each eye must not exceed the allowable tolerance for Horizontal Prism listed in the Quality Standards on the EIP job aid.</p> <p>2. Horizontal Prism: The sum of the horizontal prism readings (In & Out) for each eye must not exceed the allowable tolerance for Horizontal Prism listed in the Quality Standards on the EIP job aid.</p>		<p>Vertical Prism</p> <p>3. Vertical Prism: The sum of the vertical prism readings (Up & Down) for each eye must not exceed the allowable tolerance for Vertical Prism listed in the Quality Standards on the EIP job aid.</p> <p>4. Vertical Prism: The sum of the vertical prism readings (Up & Down) for each eye must not exceed the allowable tolerance for Vertical Prism listed in the Quality Standards on the EIP job aid.</p>	
<p>Vertical Imbalance</p> <p>5. Vertical Imbalance: The sum of the vertical prism readings (Up & Down) for each eye must not exceed the allowable tolerance for Vertical Imbalance listed in the Quality Standards on the EIP job aid.</p> <p>6. Vertical Imbalance: The sum of the vertical prism readings (Up & Down) for each eye must not exceed the allowable tolerance for Vertical Imbalance listed in the Quality Standards on the EIP job aid.</p>		<p>Horizontal Imbalance</p> <p>7. Horizontal Imbalance: The sum of the horizontal prism readings (In & Out) for each eye must not exceed the allowable tolerance for Horizontal Imbalance listed in the Quality Standards on the EIP job aid.</p> <p>8. Horizontal Imbalance: The sum of the horizontal prism readings (In & Out) for each eye must not exceed the allowable tolerance for Horizontal Imbalance listed in the Quality Standards on the EIP job aid.</p>	





Note: Use the "Inspection Tolerance" for optimized prescriptions from RxO.

INSP. TOLERANCE							
Sphere	Tol	Cyl	Tol	Axis	Tol	Prism 1	Prism 2
R: +1.75	0.13	-2.00	0.13	180	2	0.12I	0.37D
L: +2.00	0.13	-2.00	0.13	170	2	0.12I	0.37D

1. Set the lens meter to the single vision setting:

1A: If needed, use the Measurement selection icon to change lens type. 


1B: Press the Clear icon to erase previous readings. 

1C: Press the Lens icon to change from single lens to complete pair.  



2. Place the right eye over the lens stop:

2A: Position the bottom of the frame so its resting against the frame table with the right-side nose pad/bridge in contact with the left nose pad that's attached to the frame table.



3. Move the cross to the center of the target:

3A: While keeping the frame against the frame table, center the cross within the target. Use your left hand to move the cross right or left and the frame table lever to move it up or down.


3B: Make fine movements to get the cross in the exact center of the target. As the cross approaches the center, the screen will display Alignment OK. When the cross is centered it will change color from pink to blue and the screen will display Marking OK.



4. Clamp the lens:


4A: After the cross is centered in the target and the screen displays Marking OK, clamp the R lens in place.

4B: When the lens is clamped, the cross may move a little out of alignment and need to be readjusted. With one hand, release just a bit of pressure from the lens clamp, and readjust the alignment cross so it's perfectly centered and re-clamp.




5. Press the Memory Button:

5A: When the lens is perfectly placed and clamped, press the Memory Button. The color of the measurement value box will change, indicating the measurement values are saved.




6. Unclamp and position the lens over the Bifocal:

6A: Release the lens clamp and use the frame table lever to move the table towards yourself until the Bifocal is centered over the lens stop, then re-clamp the lens.




7. Press the Memory button to capture the Add power:

7A: Make sure the lens is clamped and the frame is resting against the table before capturing the Add power, unclamped lenses will read inaccurately.




8. Remove the right lens:

8A: Lift the lens clamp and remove the R lens.



INSP . TOLERANCE				
Sphere	Tol	Cyl	Tol	Axis
R: + 1.75	0.13	- 2.00	0.13	180
L: + 2.00	0.13	- 2.00	0.13	170

[illegible]

	Power (@	
Id	90	180
	- 4.42	- 4.33
	- 4.99	- 4.51


Vendor

11-1a. Checking Lenses over 2.75D

11-1A: Compare the value listed on the Lab Ticket to the value you read using the EIP Alignment Grid:
 11-1B: Is the Right OC within 1.5 mm of the prescribed value? If YES, the lens passes.
 11-1C: Is the Left lens within 1.5 mm of the prescribed value? If YES, the lens passes.
 11-1D: Is the pair within a total of 2.5 mm from prescribed? If YES, the eyewear passes inspection for PD. Continue to the Check OC Height.

11-1b. Checking Distance PD in lenses under 2.75 D:

11-1A: 1. Use an AR pen to mark where the OC should be both horizontally and vertically.
 For example- the lens below marked in the lensometer at 30 mm. The lab ticket listed 32 as the correct PD. Mark a cross at 32 mm




11-1B: Center the cross in the Lens Stop. Use the Axis Marking Pins to ensure the middle pin is directly over the center of the cross.
 11-1C: Clamp the lens and press the memory button to capture the prism value.
 11-1D: The horizontal prism values may be up to 0.33D in each eye.
 11-1E: If both lenses are equal to or less than 0.33, for a total of 0.67, the pair passes inspection for PD



In tolerance
below .33D

12. The OC's were spotted in Step 1. If needed re-spot by:

- 12A: Align the Right Lens until the Target is centered and the screen reads "Marking Ok"
- 12B: Use the marking pins to mark the Right OC
- 12C: Align the Left lens until the Target is centered and the screen reads "Marking OK"
- 12D: Use the marking pins to mark the Left OC



13. Measure the OC Height Using the EIP JobAid:



13A: Use the Alignment Grid on the EIP Job aid by aligning the frame to the dark bottom line and reading off the scale.

Note: the image to the right is an example of a Right eye with the red OC mark reading 20 mm on the scale.

13B: Determine the values on the Right and Left Lenses.

13C: Compare to the work ticket- If the values match the pair passes.


If the values do not match continue to Near PD.



13-1. Determine the power of the lenses at 90°:

13-1A: Use the Lab Ticket to determine the power of the lens at 90- The image on the right shows the Right lens is over 3.75. Go to the next step.

13-1B: If the power is less than 3.75 go to Checking Lenses under 3.75




13-1a. Checking Lenses over 3.75D :

13-1A: Compare the value listed on the Lab Ticket to the value you read using the EIP Alignment Grid:

13-1B: Is the Right OC within 1.0 mm of the prescribed value? If YES, the lens passes.

13-1C: Is the Left lens within 1.0 mm? If YES, the lens passes.


13-1D: Are they within 1.0 mm of each other? If YES, the pair passes inspection for OC Height



13-1b. Checking OC's in lenses under 3.75 D:

13-1A: Use an AR pen to mark where the OC should be both horizontally and vertically.

For example- the lens below marked in the lensometer at 18 mm. The lab ticket listed 20 as the correct OC. Mark a cross at 20 mm.




13-1B: Center the cross in the Lens Meter. Use the Axis Marking pins to ensure the middle pin is directly over the center of the cross.

13-1C: Clamp the lens and press the memory button to capture the prism value.

13-1D: The values may be up to 0.33D in each eye.

13-1E: If both lenses have prism check the direction. If both down or both up and each lens passed individually the pair passes. If one lens has UP prism and one DOWN- add the values together, the result must be equal to or less than 0.33D.



14. Measure the Near PD Using the EIP JobAid:


14A: Use the EIP job aid to center the frame on the grid with the outer left edge of the Left Lens seg on the zero line.

Note: The image to the right is an example of a Right eye with the white fitting cross reading 60 mm on the scale.

14B: Read the measurement on the inner most left edge of the Right lens for the near PD.

14C: Compare values to the work ticket- If the values match, continue to checking Seg Height.

If the values do not match continue to the next step



LAB		Printed: 1/19/2018 1:03:31 PM	Customer: ROMERO, FAUSTO	Associate: HERNAN, HENRY	Station: 1633	Order #: 858383609	Text
Lens:		Frame: 44287800013 - Overlight 50 Str Pgm CGS STR PGM L10					
Material: Plastic Bifocal Straight Top 5728 0004 13		Source: Manual		Type: Normal		Total: 48.50 32.0 36.0 64.0 180.0 18.00	
Rx:		Sphere		Col	Axis	Axis	Axis
R: +1.75		-0.75	85	0.00	0.00	0.00	0.00
L: +1.75		0.00	0	0.00	0.00	0.00	0.00
Cyl:		Axis		Axis		Axis	
R: 0.00		0.00		0.00		0.00	
L: 0.00		0.00		0.00		0.00	
Add:		Axis		Axis		Axis	
R: 0.00		0.00		0.00		0.00	
L: 0.00		0.00		0.00		0.00	
Seg Ht:		OC-H		OC-H		OC-H	
R: 11.00		16.00		13.0		+1.75 +2.48	
L: 11.00		16.00		13.0		+1.75 +1.75	
Seg Ht:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		1 Hour	
L: 4.50		-4.00		-4.00		Edger Setting	
Lens Add Dia:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Lab Inspector	
L: 4.50		-4.00		-4.00		Retail Inspector	
Time Out:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Time Out	
L: 4.50		-4.00		-4.00		Time Out	
Carder ID: 3942		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Time Out	
L: 4.50		-4.00		-4.00		Time Out	

14-1. Determine the allowable Near PD tolerance.

14-1A: Compare the value listed on the Lab Ticket to the value you read using the EIP JobAid:

14-1B: Is the Right lens within 1.5mm of specified PD? If YES, the lens passes.

14-1C: Is the Left lens within 1.5 mm? If YES, the lens passes.

14-1D: Is the pair within a total of 2.5 mm from prescribed? If YES, the eyewear passes inspection for PD. Continue to the Check Seg Height.

LAB		Printed: 1/19/2018 1:03:31 PM	Customer: ROMERO, FAUSTO	Associate: HERNAN, HENRY	Station: 1633	Order #: 858383609	Text
Lens:		Frame: 44287800013 - Overlight 50 Str Pgm CGS STR PGM L10					
Material: Plastic Bifocal Straight Top 5728 0004 13		Source: Manual		Type: Normal		Total: 48.50 32.0 36.0 64.0 180.0 18.00	
Rx:		Sphere		Col	Axis	Axis	Axis
R: +1.75		-0.75	85	0.00	0.00	0.00	0.00
L: +1.75		0.00	0	0.00	0.00	0.00	0.00
Cyl:		Axis		Axis		Axis	
R: 0.00		0.00		0.00		0.00	
L: 0.00		0.00		0.00		0.00	
Add:		Axis		Axis		Axis	
R: 0.00		0.00		0.00		0.00	
L: 0.00		0.00		0.00		0.00	
Seg Ht:		OC-H		OC-H		OC-H	
R: 11.00		16.00		13.0		+1.75 +2.48	
L: 11.00		16.00		13.0		+1.75 +1.75	
Seg Ht:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		1 Hour	
L: 4.50		-4.00		-4.00		Edger Setting	
Lens Add Dia:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Lab Inspector	
L: 4.50		-4.00		-4.00		Retail Inspector	
Time Out:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Time Out	
L: 4.50		-4.00		-4.00		Time Out	
Carder ID: 3942		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Time Out	
L: 4.50		-4.00		-4.00		Time Out	

15. Measure the Seg Height Using the EIP JobAid:


15A: Use the Alignment Grid on the EIP Job aid by aligning the frame to the dark bottom line and reading off the scale.

Note: the image to the right is an example of a Right eye with the Seg Ht. mark reading 16.5 mm on the scale.

15B: Determine the values on the Right and Left Lenses.

15C: Compare to the work ticket- If the values match the pair passes.

If the values do not match continue to the next step



LAB		Printed: 1/19/2018 1:03:31 PM	Customer: ROMERO, FAUSTO	Associate: HERNAN, HENRY	Station: 1633	Order #: 858383609	Text
Lens:		Frame: 44287800013 - Overlight 50 Str Pgm CGS STR PGM L10					
Material: Plastic Bifocal Straight Top 5728 0004 13		Source: Manual		Type: Normal		Total: 48.50 32.0 36.0 64.0 180.0 18.00	
Rx:		Sphere		Col	Axis	Axis	Axis
R: +1.75		-0.75	85	0.00	0.00	0.00	0.00
L: +1.75		0.00	0	0.00	0.00	0.00	0.00
Cyl:		Axis		Axis		Axis	
R: 0.00		0.00		0.00		0.00	
L: 0.00		0.00		0.00		0.00	
Add:		Axis		Axis		Axis	
R: 0.00		0.00		0.00		0.00	
L: 0.00		0.00		0.00		0.00	
Seg Ht:		OC-H		OC-H		OC-H	
R: 11.00		16.00		13.0		+1.75 +2.48	
L: 11.00		16.00		13.0		+1.75 +1.75	
Seg Ht:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		1 Hour	
L: 4.50		-4.00		-4.00		Edger Setting	
Lens Add Dia:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Lab Inspector	
L: 4.50		-4.00		-4.00		Retail Inspector	
Time Out:		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Time Out	
L: 4.50		-4.00		-4.00		Time Out	
Carder ID: 3942		Axis		Axis		Axis	
R: 4.50		-4.00		-4.00		Time Out	
L: 4.50		-4.00		-4.00		Time Out	

INSP. TOLERANCE							
Sphere	Tol	Cyl	Tol	Axis	Tol	Prism 1	Prism 2
R: +1.75	0.13	-2.00	0.13	180	2	0.12I 0.37D	
L: +2.00	0.13	-2.00	0.13	170	2	0.12I 0.37D	