

The new invisible hearing solution

nuance audio
in tune with life

OUR MANIFESTO


Empowering
the world to hear,
see, and connect

Tune into a world where the beauty of sound knows no limits, where you can hear every detail and immerse yourself in the richness of life.

Because the better you hear, the more you can connect with what matters.

CONTEXT

~73M Total Americans with hearing loss



Source: Global Burden of Disease Study (GBD; Haile et. al., 2024)

Key barriers to hearing aid adoption:

Stigma

Comfort

Price

Source: Reasons for the non-use of hearing aids: perspectives of non-users, past users, and family members, Franks & Timmer, 2024

BREAKTHROUGH INNOVATION

Hear better in
noisy environments
with the all-in-one
hearing glasses

The new invisible open-ear OTC hearing solution for perceived mild to moderate hearing loss.

Hear what you see
with directional
microphones

Hear comfortably
with discreet
open-ear speakers

What Makes Nuance Audio Hearing Glasses Unique

ALL-IN-ONE SOLUTION

A comprehensive solution addressing both vision and hearing needs in one elegant device.

INVISIBLE HEARING AID

Break the stigma with discreet, stylish glasses that seamlessly integrate vision and hearing support.

OPEN-EAR COMFORT

Lightweight, open-ear speakers offer all-day comfort and avoid pain or irritation often caused by traditional hearing aids.

DIRECTIONAL AMPLIFICATION

Enhance speech clarity in noisy environments with directional microphones & beamforming.

EASY TO USE

Simple to set up on your own with a user-friendly experience that does not require audiologist assistance.

STYLES & PACKAGING

Available in two styles, Square 54/56 and Panthos, and two different colors to complement any look.

Inside the packaging you will find:

- 1 x Nuance Audio Glasses
- 1 x Glasses Case
- 1 x Cleaning Cloth
- 1 x Charging Pad, with cable included



ACCESSORIES


Designed to be discreet and easy to use while maintaining the same sophisticated look as the eyewear.



LENSES

Nuance Audio Glasses are best paired with Transitions® lenses for seamless hearing and sight from day to night.

Available in non-prescription lenses in Grey and Emerald at purchase, and in a wide color range for prescription lenses on demand.



Option to add your preferred optical lenses based on your vision needs.

Transitions and the Transitions logo are registered trademarks of Transitions Optical, Inc. used under license by Transitions Optical Limited. Photochromic performance and polarization is influenced by temperature, UV exposure and lens material.

USER EXPERIENCE


Nuance Audio App: The essential Nuance Audio companion

 Download for iOS

 Download for Android

Personalize the Nuance Audio experience with the Nuance Audio app. With the app users can:

- 1 Manage noise reduction level
- 2 Switch between alternative audio modes: Frontal vs. All-Around
- 3 Set your preferred preset options for a tailored hearing experience



The Nuance Audio App is compatible with iOS devices running iOS 16 or later (release 0.3.7 onward) and Android devices running Android 12 or later with SDK 24 (release 0.4.1 onward).