

INSTRUCTIONS FOR CONVERTING PORTABLE TO SIMPLIFIED DEMO SYSTEM

Step 1

Open the Nuance Audio portable solution



Step 2

Unplug and remove the neckband



Step 3

Unplug and remove the smartphone



Step 4

Unplug and remove the Nuance Audio charging pad



Step 5

Now the tray is empty. Please lift and remove it



Step 6

Unplug cables A, B, D from the powerbank




Step 7

Remove and un-clip cable C following its path. Pull that out from the powerhub and peel off any glue residue



Step 8

Remove and put aside cable A



FINAL CONFIGURATION

Here is the final selection of items to use in order to perform the Nuance Audio demo experience

